


HD CSPS SYSTEM

CONTAINER / STERILE PACK HIGH DENSITY STORAGE SYSTEM



MAIN FEATURES:

- Stainless Steel container/pack beds will not rust and are scratch resistant. Unlike plastic or powder coated shelves, our stainless steel beds will not release unwanted particles into your sterile environment.
- A complete modular system: spacing between shelves can be reconfigured by the end user, without using any tools
- Access to entire Container/Pack via 110% retractable extension slides
- Beds have raised side and back walls to assure containers/trays will not fall while beds are moved in and out
- Load capacity per shelf: 37kg/88lbs



SS-U-CNT-SHELF-HD

10 units per box, telescopic slides and handles should be ordered separately

The container beds can adjusted every 1"


CONTACT US FOR MORE DETAILS

+1 888-276-4750 | info@pegasusmedical.net | www.pegasusmedical.net

RACK DIMENSIONS




COLOR CODED ERGONOMIC HANDLE SYSTEM



Matte aluminum # HNDL-112MM-ALU Red # HNDL-112MM-RED
Blue # HNDL-112MM-BLUE Yellow # HNDL-112MM-YLW
(10 pcs handles per pack)

U Type



* Racks can be ordered with leveling feet instead of wheels

CONTACT US FOR MORE DETAILS

V200119

+1 888-276-4750 | info@pegasusmedical.net | www.pegasusmedical.net