

APOLLO CARTS

THE CHAMELEON SOLUTION. ADAPTING
TO EVERY HEALTHCARE CHALLENGE



SUTURES

CATH LAB

GENERAL SUPPLY

MEDICATION

STERILE CORE

APOLLO CARTS

THE FUTURE OF HEALTHCARE

Healthcare is evolving, becoming leaner, more efficient, and cost-effective. Yet, bulky, outdated cabinets keep hospitals stuck in rigid setups. The following are some of the challenges healthcare facilities are facing as they adapt to modern demands.



CHANGING DEMOGRAPHICS

Hospitals are adapting spaces to reflect the cultural and health needs of their local communities



TECH ADVANCEMENTS

Rooms must stay flexible to accommodate constantly evolving medical technology



OUTPATIENT SERVICES

Care is shifting toward faster, more accessible outpatient treatments



INFECTION CONTROL

Designs now prioritize hygiene with antimicrobial surfaces and easy to clean



DOWNSIZING HOSPITAL

Smaller, smarter facilities focus on efficiency and multi-use spaces




Apollo Storage Carts offer modular mobility and impermeable materials for a cleaner, adaptable future. Ditch the static, move forward with Apollo



SAY GOODBYE TO BULKY CABINETS

APOLLO CARTS THE FUTURE OF STORAGE!



Scratch = Bacteria & Mold

⊗ LIMITED ACCESS

Fixed cabinets restrict access and require more physical effort to reach supplies, increasing strain on nurses and slowing down workflows

⊗ RIGID CONSTRAINTS


When hospital administrations change, fixed storage solutions often become obstacles, as cabinets and layouts are locked into designs that can't easily be reassigned or rearranged to fit new priorities

⊗ FIXED & INFLEXIBLE

Permanently attached cabinets restrict room layout changes, making it difficult to adapt to evolving workflows and complicating room cleaning

⊗ ABSORBENT

Wood is rigid and prone to cracking, increasing moisture and bacteria absorption, creating a breeding ground for contamination



Scratches don't break impermeable surface

✓ LIGHTWEIGHT & NURSE-FRIENDLY

Ergonomic and mobility reduce physical strain, and allows nurses to quickly reposition carts and access supplies with ease

✓ STANDARD CORE

Apollo's modular design, based on standardized components, allows for seamless adjustments, making it easy for new administrations to quickly adapt storage to meet evolving healthcare needs

✓ FLEXIBLE & MOVABLE

Easily repositioned or reassigned, Apollo carts adapt to evolving hospital layouts and workflows, while their mobility allows for thorough room cleaning

✓ NON ABSORBENT

Polymer and aluminum repel water and microbes, offering impact-resistant flexibility and lasting integrity for a cleaner, safer environment

10 REASONS WHY APOLLO IS THE FUTURE OF STORAGE

TAMBOUR DOOR

A back-panel system lets the door slide behind modules without interference, enabling a compact, counterbalanced design without adding height or losing internal space

POLYMER SKIN

A tough, impact-resistant exterior that's easy to clean and corrosion-resistant

MOBILITY

5" casters for easy, smooth movement and repositioning

MODULAR

A flexible design that allows easy customization to meet changing needs

ALUMINUM SKELETON

A lightweight, durable aluminum frame that ensures strength and stability

DRAIN FLOOR

A floor design with drainage to prevent fluid buildup and maintain hygiene



4

DOORS

4

LOCKS

5

HEIGHTS

6

WIDTHS

THE FUTURE IS NOW

CREATE STORAGE YOUR WAY!

Discover a smarter way to store, organize, and adapt. Explore how Apollo has solutions designed for every hospital need



GENERAL SUPPLY

Versatile storage for essential hospital supplies

Page 6



CATH LAB

Organized, adaptable catheter storage

Page 8



MEDICATION

Organized for efficient med management

Page 10



SUTURES

Dedicated modules for easy identification and access

Page 12



STERILE CORE

Dedicated SPD/CSSD storage

Page 14



ANTIMICROBIAL

Extra protection against bacteria


Page 16



MRI

Non-magnetic, safe for MRI zones

Page 17



UNDER THE COUNTER

Compact storage for tight spaces

Coming Soon

APOLLO GENERAL SUPPLY

BUILT TO ADAPT, READY FOR CHANGE



ENHANCED WORKFLOW

Adjustable configurations simplify supply storage and streamline daily tasks, saving time and effort



STANDARDIZED SOLUTION

HTM71 (ISO 6040) standard ensures seamless integration with hospital storage systems, improving workflow and organization



VERSATILE SYSTEMS

Modular design adapts to any department, offering multi-purpose storage that can be reassigned as healthcare demands evolve

Designed for any department, Apollo carts follow the ISO 6040 standard for seamless integration. Their modular design allows easy customization and smooth transfers, improving workflow and efficiency



DIVIDER SYSTEM

Improve inventory rotation and reduce waste by organizing supplies with our Kanban work method



LABELING SYSTEM

Reduce errors and improve efficiency by adding identifying labels to the storage structure



COLOR CODING SYSTEM

Cut supply retrieval time by up to 25% with a visual system that organizes and identifies items at a glance



DEDICATED MODULES

Dedicated solutions to all your storage needs, including antibacterial, clear, CSSD/SPD modules, and more

PRECONFIGURED APOLLO GENERAL SUPPLY



APL-1 NARROW 72.4\"
GENERAL SUPPLY CART

- 1x Shelves
- 2x 5 cm (2\") Tray
- 5x 10 cm (4\") Tray
- 2x 20 cm (8\") Tray
- 9x Color-Coding Handle & Strip
- 18x Ver. Label Holders & Cover
- 2x Vertical Push Handles
- 1x VF Dust Cover



H - 184 cm | 72.4\"
D - 73 cm | 28.6\"
W - 55 cm | 21.6\"



APL-1 WIDE 72.4\"
GENERAL SUPPLY CART

- 2x 5 cm (2\") Basket
- 6x 10 cm (4\") Basket
- 2x 20 cm (8\") Basket
- 9x Color-Coding Handle & Strip
- 2x Vertical Push Handles
- 1x VF Dust Cover



H - 184 cm | 72.4\"
D - 53 cm | 20.7\"
W - 75 cm | 29.5\"



APL-2 NARROW 72.4\"
GENERAL SUPPLY CART

- 4x Shelves
- 2x 5 cm (2\") Tray
- 10x 10 cm (4\") Tray
- 2x 20 cm (8\") Tray
- 14x Color-Coding Handle & Strip
- 28x Ver. Label Holders & Cover
- 2x Vertical Push Handles
- 1x VF Dust Cover



H - 184 cm | 72.4\"
D - 73 cm | 28.6\"
W - 106 cm | 41.7\"



APL-2 WIDE 72.4\"
GENERAL SUPPLY CART

- 8x Shelves
- 2x 5 cm (2\") Tray
- 4x 10 cm (4\") Tray
- 2x 20 cm (8\") Tray
- 8x Color-Coding Handle & Strip
- 2x Vertical Push Handles



H - 184 cm | 72.4\"
D - 53 cm | 20.7\"
W - 146 cm | 57.6\"



APL-3 NARROW 72.4\"
LINNEN & SUPPLY CART

- 14x Shelves
- 2x Vertical Push Handles




H - 184 cm | 72.4\"
D - 73 cm | 28.6\"
W - 157 cm | 61.7\"



APL-3 WIDE 72.4\"
GENERAL SUPPLY CART

- 6x Shelves
- 9x 10 cm (4\") Basket
- 6x 20 cm (8\") Basket
- 15x Color-Coding Handle & Strip



H - 184 cm | 72.4\"
D - 53 cm | 20.7\"
W - 217 cm | 85.4\"

APOLLO CATH LAB

STREAMLINING YOUR WORKFLOW



Apollo Cath Lab Carts accommodate all catheter sizes, keeping supplies organized, visible, and accessible. Customize with hooks, baskets, and adjustable shelves for a seamless workflow



Catheter security, without compromise



4 locking systems to ensure controlled access and maintain compliance



FIT FOR ALL HEIGHTS

Our modular design easily adapts to the varying heights of catheters, optimizing storage efficiency



DUAL STORAGE

Catheters and guides are stored in pouches and boxes, requiring different ways to keep them organized. Use our dedicated hooks and baskets for easy access



ALL-IN-ONE PLACE

Keep all cath lab supplies in a single cart with a mix & match of modules, ensuring quick and easy access



PRECONFIGURED APOLLO CATH LAB



APL-1 NARROW 72.4" CATH & SUPPLY CART

- 1x Shelf
- 3x 9 Hook Cath Pull-outs
- 2x 10 cm (4") Tray
- 1x 20 cm (8") Tray
- 3x Color-Coding Handle & Strip
- 6x Ver. Label Holders & Cover
- 5x Side Cath Box
- 2x Side Acc. Rails
- 2x Vertical Push Handles



H - 184 cm | 72.4"
D - 73 cm | 28.6"
W - 55 cm | 21.6"



APL-1 WIDE 72.4" 10 CATH CART

- 2x Shelves
- 10x 6 Hook Cath Pull-outs
- 1x VF Dust Cover



H - 184 cm | 72.4"
D - 53 cm | 20.7"
W - 75 cm | 29.5"



APL-2 NARROW 72.4" BOXED CATH CART

- 2x Shelves
- 2x Boxed Cath Pull-outs
- 2x Hor. Label Holders & Cover
- 2x Vertical Push Handles
- 1x VF Dust Cover



H - 184 cm | 72.4"
D - 73 cm | 28.6"
W - 106 cm | 41.7"



APL-3 NARROW 72.4" 6 CATH & SUPPLY CART

- 3x Shelves
- 3x 9 Hook Cath Pull-outs
- 6x 10 cm (4") Tray
- 4x 20 cm (8") Tray
- 10x Color-Coding Handle & Strip
- 20x Ver. Label Holders & Cover
- 2x Vertical Push Handles



H - 184 cm | 72.4"
D - 73 cm | 28.6"
W - 157 cm | 61.7"



APL-2 NARROW 72.4" ENDOVASCULAR CART

- 3x Shelves
- 3x 9 Hook Cath Pull-outs
- 1x 5 cm (2") Tray
- 4x 10 cm (4") Tray
- 3x 20 cm (8") Tray
- 8x Color-Coding Handle & Strip
- 16x Ver. Label Holders & Cover
- 2x Vertical Push Handles
- 1x VF Dust Cover



H - 184 cm | 72.4"
D - 73 cm | 28.6"
W - 106 cm | 41.7"



APL-3 NARROW 72.4" CYSTO CART

- 8x Shelves
- 3x 9 Hook Cath Pull-outs
- 6x 10 cm (4") Tray
- 4x 20 cm (8") Tray
- 10x Color-Coding Handle & Strip
- 20x Ver. Label Holders & Cover
- 2x Vertical Push Handles



H - 184 cm | 72.4"
D - 73 cm | 28.6"
W - 157 cm | 61.7"

APOLLO MEDICATION

SIMPLIFIED & EFFICIENT STORAGE



Apollo Medication Carts enhance workflow with two storage options: a narrow, concealed unit for secure, discreet transport between departments, and a wide, open-access unit for quick, convenient use in clinical workspaces, pharmacies, and medication rooms



SEAMLESS FLOW

ISO 6040 standardized cassettes easily move between racks, Apollo carts, and procedure carts, streamlining workflow



SMART TRANSPORT

Narrow front keep medication info private and easy to navigate in tight spaces, while wide front offer full visibility for high-turnover meds



SECURE MEDICATION HANDLING

Reduce risks of misuse and errors by ensuring safe and compliant handling of medications, with customizable medication cassettes



3 level cassette = 1 panel

For ease of use and interchangeability, we recommend the 3 level Medicine Cassette

PRECONFIGURED APOLLO MEDICATION



APL-1 NARROW 72.4" MED TRANSPORT

4x 3 Level MedCast
2x Vertical Push Handles
1x VF Dust Cover



H - 184 cm | 72.4"
D - 73 cm | 28.6"
W - 55 cm | 21.6"



APL-3 NARROW 72.4" MED PRIVACY

12x 3 Level MedCast
2x Vertical Push Handles



H - 184 cm | 72.4"
D - 73 cm | 28.6"
W - 157 cm | 61.7"



APL-1 WIDE 72.4" MULTIPURPOSE MED

2x 3 Level MedCast
1x Drawer Front on Telescopic Rails - TDF-C2
1x 5 cm (3") Tray
2x 10 cm (4") Tray
2x 20 cm (9") Tray
1x VF Dust Cover



H - 184 cm | 72.4"
D - 53 cm | 20.7"
W - 75 cm | 29.5"



APL-2 WIDE 72.4" MED STORAGE

8x 3 Level MedCast
2x Vertical Push Handles



H - 184 cm | 72.4"
D - 53 cm | 20.7"
W - 146 cm | 57.6"

The empty cassettes return to the pharmacy



APL-2 WIDE 43" MED STORAGE

4x 3 Level MedCast
2x Key Locks



H - 109 cm | 43"
D - 53 cm | 20.7"
W - 146 cm | 57.6"



The new units of medication cassettes arrive filled from the central pharmacy and are distributed to the wards where the cassettes have been exchanged (empty for full)

APOLLO SUTURES

EFFORTLESS ORGANIZATION AND PROTECTION



Designed for surgical precision, Apollo Suture carts offer two storage options: pull-out shelves for quick access and transparent trays for organized color-coded backup. The ideal setup combines both, with shelves on top and trays below for efficient use in one cart



**Protect Every Suture,
Prevent Every Loss**

4 locking systems to reduce waste
and prevent unauthorized access



DEDICATED STORAGE

Telescopic side shelf fits 4-8 rows of suture boxes for quick access



CRYSTAL CLEAR COLOR CODING

Clear trays allows easy organization of color-coded suture boxes for fast identification



PROTECTED & VISIBLE

Glass doors offer visibility while keeping suture boxes protected from dust and unauthorized access when used with electronic locks




PRECONFIGURED APOLLO SUTURES



APL-2 NARROW 72.4\"/>

2x Pull Out Suture Shelves Set
2x Suture Shelf Kit
10x 5 cm (3\") Clear Tray
2x Vertical Push Handles
1x VF Dust Cover




H - 184 cm | 72.4\"/>



APL-2 WIDE 72.4\"/>

22x 5 cm (3\") Clear Tray
2x Vertical Push Handles




H - 184 cm | 72.4\"/>



APL-3 NARROW 72.4\"/>

2x Pull Out Suture Shelves Set
2x Suture Shelf Kit
22x 5 cm (3\") Clear Tray
2x Vertical Push Handles




H - 184 cm | 72.4\"/>



APL-1 NARROW 72.4\"/>

1x Pull Out Suture Shelves Set
1x Suture Shelf Kit
5x 5 cm (3\") Clear Tray
2x Vertical Push Handles
1x VF Dust Cover




H - 184 cm | 72.4\"/>



APL-2 NARROW 72.4\"/>

12x HPL Shelves
4x 20 cm (8\") Tray
4x Color-Coding Handle & Strip
8x Ver. Label Holders & Cover
2x Vertical Push Handles
1x VF Dust Cover




H - 184 cm | 72.4\"/>



APL-1 WIDE 72.4\"/>

11x 5 cm (3\") Clear Tray
1x VF Dust Cover



H - 184 cm | 72.4\"/>

APOLLO STERILE CORE

SMART STORAGE & SAFE TRANSPORT



CLOSED-LOOP DISTRIBUTION

One cart does it all: store in the sterile area, transport to the procedure room, and keep it nearby during use



PROTECTED

Fully enclosed design ensures safe transport and prevents cross-contamination in single-corridor setups



NO-TOUCH METHOD

Stainless steel modules support no-touch handling, while telescopic slides manage containers with heavy loads with ease



PAY ATTENTION! DEDICATED MODULES FOR STERILE CORE


Please follow our best practice guidelines for using our sterile core modules to keep supplies organized, protected and accessible while preventing contamination

PRECONFIGURED APOLLO STERILE CORE




APL-2 NARROW 72.4" FOR NO TOUCH & CONTAINERS

8x No Touch Wire Shelf
8x Container Wire Shelf
16x Non-Slip Silicone Mat
2x Vertical Push Handles
1x VF Dust Cover




APL-2 NARROW SPD WORKSTATION

8x No Touch Wire Shelf
8x Non-Slip Silicone Mat
6x 10 cm (4") Half Tray - 30x40
6x Color-Coding Handle & Strip




APL-1 NARROW 72.4" SPD TRANSPORT

4x No Touch Wire Shelf
2x Container Wire Shelf
2x Container Solid Shelf
2x Solid Shelf Handle
8x Non-Slip Silicone Mat
2x Vertical Push Handles
1x VF Dust Cover




APL-1 WIDE 72.4" MODULAR SHELF SPD CART

9x HPL Shelves
2x Vertical Push Handles
1x VF Dust Cover




APL-2 NARROW 43" FOR CONTAINERS

8x Container Wire Shelf
8x Non-Slip Silicone Mat
1x Horizontal Push Handles
1x VF Dust Cover



APL-2 NARROW 57.6" FOR CONTAINERS

6x Container Wire Shelf
6x Container Solid Shelf
6x Red Solid Shelf Handle
12x Non-Slip Silicone Mat
2x Vertical Push Handles
1x VF Dust Cover



DEFEND AGAINST INFECTIONS WITH APOLLO + ANTIMICROBIAL TRAYS



24H ANTIMICROBIAL PROTECTION

Storage modules are made with antimicrobial materials to help reduce contamination risks



FLEXIBLE TRAY OPTIONS

Available in four tray sizes to fit your storage needs and keep everything organized



CLEANING STANDARD COMPLIANT

Antimicrobial materials and durable surfaces withstand harsh hospital-grade cleaners for long-lasting hygiene



Preventing secondary infections in procedure rooms and ORs starts with smarter storage. Our carts reduce contamination risks with antimicrobial modules, smooth surfaces for easy cleaning, and secure storage, ensuring a compliant environment for safer patient care

DON'T DISRUPT THE MAGNETIC FIELD

APOLLO MRI CERTIFIED



MR CONDITIONAL CERTIFIED

Apollo is tested and approved to ensure it won't interfere with MRI machines or pose any safety risks when used under specified conditions



NON-MAGNETIC

Crafted from non-magnetic metals, our products ensure safety and reliability within the MRI environment, without any magnetic interference

MRI Safety Information



MR Conditional

Non-clinical testing has demonstrated the "Apollo cart" is MR Conditional. The device may be safely used in the MR environment under the following conditions. Failure to follow these conditions may result in injury to the patient and/or item malfunction.

| Parameter | Condition of Use/Information |
|--|----------------------------------|
| Static Magnetic Field Strength (B_0) | 1.5 T, 3 T |
| Maximum Spatial Field Gradient | 0.2 T/m (20 gauss/cm) |
| Magnetic Field Safety Threshold | Do not exceed 60 mT (600 gauss). |

If information about a specific parameter is not included, there are no conditions associated with that parameter.



DESIGN YOUR IDEAL CART

STEP 1 - MODEL

WIDE (E) TYPE



VISIBILITY & ACCESS

Wide front design makes supplies easy to see and reach



SPACE-SAVING DEPTH

Shallow profile preserves corridor space and improves traffic flow



HIGH-VOLUME READY

Designed to handle frequent use and high turnover with easy loading and retrieval



NARROW (U) TYPE



COMPACT & MANEUVERABLE

Slim profile makes it easy to move and navigate in tight spaces



SPACE-SAVING CAPACITY

Small footprint increases storage while saving floor space





EFFICIENT STORAGE

Maximizes capacity, making it a great cabinet replacement for limited spaces





STEP 2 - WIDTH



SINGLE



DOUBLE



TRIPLE



STEP 3 - HEIGHT



209 cm | 82.3" – APL-4.6

184 cm | 72.4" – APL-4.0

146 cm | 57.6" – APL-3.0

109 cm | 43" – APL-2.0

STEP 4 - CASTER



5" Wheel

Smooth mobility for effortless maneuvering



4" Wheel

2 cm shorter, similar mobility



Level Caster

Stable like a cabinet, movable for cleaning (-7.5 cm height)



Level Feet

Fixed cabinet setup (-13 cm height)

STEP 5 - SLIDERS


Tilt Panels

Modules are stored horizontally and can be tilted while pulled out




Telescopic Slides

Designed to securely handle heavy loads with full extension for easy access



You can use a combination of Tilt Panels and Telescopic Slides in the following %:



STEP 6 - DOORS



Tempered Glass (GD)

Easy access with clear visibility



Tambour (TD)

Small footprint with no swing clearance



Solid (SD)

Secure storage with full privacy



No Door (ND)

Open access for high-turnover supplies

STEP 7 - LOCKING SYSTEM



Key Lock

Basic security
APL-ALD-KEYLCK



Electronic

No keys to lose
APL-ALD-CODELCK (GD, SD)
APL-TD-CODELCK (TD)



Local Management RIFD

ID-based access
PLY-CRDCLK-LM



Central Management RIFD

User login tracking & access control
PLY-CRDCLK-CM

OPTIONAL - DRAWER CONFIGURATION (FOR E CART ONLY)


Improve workflow efficiency with easy access to frequently used items while keeping less frequently needed supplies behind doors.

Available with Glass, Solid, and Tambour doors


Height 109 cm - (2-panels)



DDF-C1 | TDF-C1



DDF-C2 | TDF-C2



DDF-C3 | TDF-C3



ACCESSORIES TO OPTIMIZE YOUR CART

ISO6040 MODULES



ABS - Cream

Trays

5 cm - T054060
10 cm - T104060
15 cm - T154060
20 cm - T204060

Baskets

10 cm - B104060
20 cm - B204060


Lid - L4060



Antimicrobial - Gray

Trays

5 cm - T054060-AM-G
10 cm - T104060-AM-G
15 cm - T154060-AM-G
20 cm - T204060-AM-G




Polycarbonate - Clear

Trays

5 cm - T054060PC
10 cm - T104060PC
20 cm - T204060PC

Baskets

10 cm - B104060PC
20 cm - B204060PC



High Resistance - Yellow

Trays

10 cm - T104060-ND-HR
20 cm - T204060-HR

Lid - L4060-HR

MAXIMUM CONTAINER CAPACITY FOR E-TYPE AND U-TYPE CARTS:

1 Panel = 9 Intervals

Panel with decimals (.)

- Decimal n° = n° of intervals

(+) Shows number of intervals left after filled with maximum allowed for each container

| | | | 2 Panels 18 Intervals | 3 Panels 27 Intervals | 4 Panels 36 Intervals | 4.6 Panels 42 Intervals |
|----------------------|--|----------------|--------------------------|--------------------------|--------------------------|----------------------------|
| 20 cm 5 Intervals | | TD GD/SD/ND | 3 (+2) 3 (+3) | 5 (+1) 5 (+2) | 7 7 (+1) | 8 (+1) 8 (+2) |
| 15 cm 4 Intervals | | TD GD/SD/ND | 4 (+1) 4 (+2) | 6 (+2) 6 (+3) | 8 (+3) 9 | 10 (+1) 10 (+2) |
| 10 cm 3 Intervals | | TD GD/SD/ND | 5 (+2) 6 | 8 (+2) 9 | 11 (+2) 12 | 13 (+2) 14 |
| 5 cm 2 Intervals | | TD GD/SD/ND | 8 (+1) 9 | 13 13 (+1) | 17 (+1) 18 | 20 (+1) 21 |


TD - Tambour door GD - Glass door SD - Solid door ND - Opened (No door)

MODULE ACCESSORIES



Handle

ABS - Cream
LBL-HDL-E-B
LBL-HDL-U-B
Antimicrobial - Gray
LBL-HDL-E-AM-G



Color-Coding Strip

Yellow - LBL-HDL-E-YP/ LBL-HDL-U-YP
Blue - LBL-HDL-E-BP/ LBL-HDL-U-BP
Green - LBL-HDL-E-GP/ LBL-HDL-U-GP
Red - LBL-HDL-E-RP/ LBL-HDL-U-RP
+ 16 more colors available in USA



Label Holders

Holder
LBL-HLD-VER
Transparent Cover
LBL-VERCOV-FL



Dividers

| Module | ABS - Cream | Antimicrobial - Gray | Polycarbonate - Clear |
|--------|-------------------------------|---|-----------------------------------|
| 5 cm | Long - D0560 Short - D0540 | Long - D0560-AM-G Short - D0540-AM-G | Long - D0560PC Short - D0540PC |
| 10 cm | Long - D1060 Short - D1040 | Long - D1060-AM-G Short - D1040-AM-G | Long - D1060PC Short - D1040PC |
| 15 cm | Long - D1560 Short - D1540 | Long - D1560-AM-G Short - D1540-AM-G | — |
| 20 cm | Long - D2060 Short - D2040 | Long - D2060-AM-G Short - D2040-AM-G | Long - D2060PC Short - D2040PC |

STAINLESS STEEL MODULES

For sterile core and O.R. storage




No Touch Wire Shelf
SS-6040WH




No Touch Wire Basket
10 cm -WS-B104060-SS-WH
20 cm -WS-B204060-SS



Container Wire Shelf
WS-INSTR-SS



Container Solid Shelf
SS-U-CNT-SHELF



Non-Slip Silicone Mat
VERS-PGS-SIMAT

MEDICINE CASSETTE



- 3 Levels** | MED-CAST-3L
2 Levels | MED-CAST-2L
4 Levels | MED-CAST-4L
5 Levels | MED-CAST-5L



Small 8cm Drawer
White - MED-DWR-80
Clear - MED-DWR-80C



Large 12cm Drawer
White - MED-DWR-120
Clear - MED-DWR-120C

SUTURE SHELVES



Pull Out Suture Shelves

SUT-SHLF-SET*
Includes:
2 Frames
8 Shelves
16 Dividers

Suture Shelf Kit


SUT-SHLF - 4 Shelves

Suture Shelf Dividers

SUT-SHLF-DIV - 8 Dividers

* REQUIRES 2 PAIRS OF TELESCOPIC SLIDES (SOLD SEPARATELY)
Telescopic slides (P Rails) and installation must be ordered with the cart.
Please specify installation location—top two panels are recommended

CATHETER MODULES



Catheter Shelves

Shelf
CATH-SHELF-E / CATH-SHELF-U
Managers
CATHH-EU / CATHH-US



Catheter Basket

CATHBX-U-TPNL
CATHBX-E-TPNL

SHELVES



Shelf with Dividers

Shelf - SHLF-POLY-H
Dividers - SHLF-DIV-PL




6mm HPL Shelf

SHLF-HPL

ACCESSORIES TO OPTIMIZE YOUR CART

DUST COVER & WORKING SURFACE

Easy-to-clean, removable top that extends your cart's lifespan



HANDLES

Ergonomic handles ensure a secure, comfortable grip for safer, easier cart maneuvering



Vertical - Recommended for carts 3 panels or taller

Wide (E) and Narrow (U) Types - APL-ALPH-U

Handle installation prep:

Right Side - APL-PH-INST-VER-R | Left Side - APL-PH-INST-VER-L

Horizontal - Recommended for carts 2 panels or shorter

Wide (E) Type - APL-ALPH-E | Narrow (U) Type - APL-ALPH-U

Handle installation prep:

Right Side - APL-PH-INST-HOR-R | Left Side - APL-PH-INST-HOR-L

COLOR-CODING SIDE PANELS

A well-designed color scheme boosts recovery, improves experience, and supports coding and navigation

| Color | Type E | Type U |
|-------------|--------------|--------------|
| Blue Poly | APL-E-SDP-B | APL-U-SDP-B |
| Red Poly | APL-E-SDP-R | APL-U-SDP-R |
| Yellow Poly | APL-E-SDP-Y | APL-U-SDP-Y |
| Pediatric | APL-E-SDP-PD | APL-U-SDP-PD |

*Special Orders can be made for different colors, adding text and back panels




SIDE ACCESSORIES

When extra storage is needed



DIN Rail - Required for side accessories

Wide (E) Type - APL-SRAIL-E
Narrow (U) Type - APL-SRAIL-U



Side Basket
SBASKET



Bottle Holder
DLBH-PP



Single Glovebox Holder
S-GBH-ABS



Triple Glovebox Holder
TR-GBH-ABS




Waste Basket
WSB-CRT



Hazard Bin Holder
HAZ-BIB-H



Side Catheter Box
S-GBH-ABS



Adjustable Oxygen Cylinder Holder
OXBH-KIT

WORKFLOW



Contact us

info@pegasusmedical.net
1-888-276-4750

1. 

Briefing

Understand and prioritize your requirements

2. 

Space Survey

Take measurements,
capture photos, and
analyze your space

3.

Tailored Solutions

Project proposals, CAD drawings,
and tailored solutions for your space

4.

Installation & Delivery:

We handle installation and delivery processes

5. 

Support:

Our representatives and we offer ongoing support for our products from purchase through delivery and beyond



**Your work
optimized,
efficient,
and safe**

NOTES:

[illegible]

APOLLO IN ACTION

REAL-WORLD CASE STUDIES ACROSS THE GLOBE

GENERAL SUPPLY



USA | MEASE COUNTRYSIDE HOSPITAL - TAMPA, FL

Installed across departments, including ORs, 50 Apollo carts support flexible workflows and quality care with their mobility and modular design



USA | UT HEALTH AUSTIN ASC - AUSTIN, TX



USA | ORLANDO HEALTH - ORLANDO, FL

A Level 1 Trauma Center uses Apollo carts in OR, General, and GI departments to support care with reliable storage and access

APOLLO IN ACTION

REAL-WORLD CASE STUDIES ACROSS THE GLOBE

CATH LAB




USA | TEXAS CHILDREN'S HOSPITAL - HOUSTON, TX

This Level 1 Pediatric Trauma Center, the nation's largest pediatric hospital, uses Apollo Cath Lab carts for streamlined access and efficient transport between prep and procedure areas



USA | TEXAS CHILDREN'S HOSPITAL - HOUSTON, TX



USA | TEXAS CHILDREN'S HOSPITAL - HOUSTON, TX

APOLLO IN ACTION

REAL-WORLD CASE STUDIES ACROSS THE GLOBE

MEDICATION



USA | NOR-LEA HOSPITAL DISTRICT - LOVINGTON, NM

Apollo carts support organized, accessible medication storage in the Nor-Lea pharmacy department, helping streamline daily workflows



USA | NOR-LEA HOSPITAL DISTRICT - LOVINGTON, NM



APOLLO IN ACTION

REAL-WORLD CASE STUDIES ACROSS THE GLOBE

STERILE CORE



ITALY | ASST LEGNANO - CUGGIONO

As part of their full SPD overhaul, West Milan health system uses Apollo carts for both transport and in-room use, combining two functions in one efficient cart



ITALY | ASST LEGNANO - CUGGIONO

